



# HBOT MASK

## - Hyperbaric Oxygen Therapy -

The Maumee Scientific respiratory face mask has been design for Hyperbaric Oxygen Therapy patients. Optimized inhale and exhale connections reduce dead space, which improves well-being of patients during HBOT sessions.



### INTENDED USE

This adult mask is to be used for spontaneously breathing patients where the application requires separate inspiratory and expiratory flows. The mask should be used in hospitals or medical environment, by individuals that have at least minimal instructions or training on the use of it as well as the device to which the mask is intended to be connected.

- **EASILY ADJUSTABLE** face mask specifically developed for Hyperbaric Oxygen Therapy usage.
- **IMPROVED BREATHING** with one way valve, which reduces dead space while the patient exhales.
- **SEPARATED INHALE AND EXHALE CONNECTIONS** for a leak free oxygen therapy.
- **REUSABLE MASK**, easy to clean, with cleaning solution, dishwasher safe.
- **THREE SIZES** to cover a wide range of patient face morphology.

### TECHNICAL SPECIFICATIONS

FACEPIECE	Thermoplastic elastomer
SET OF STRAPS WITH HARNESS	Polyester / Elastic in braided cotton
INHALE VALVE	Natural rubber
WEIGHT	117 g
PART NUMBER	MS.80.0001

CE pending.

